

Kinetic and label-free, live content imaging assays for neurite outgrowth in primary, iPSC-derived and immortalised neurons

S. L. Alcantara, P. Garay, L. McGillicuddy, T. O'Callaghan, T. Dale, D. Appledorn, V. Groppi, & D.J. Trezise
Essen BioScience, Welwyn Garden City, UK & Ann Arbor, Michigan, 48108 USA



www.essenbioscience.com

Summary & Impact

- The study of neurite dynamics is fundamental to the investigation of neuropathological disorders, neuronal injury, regeneration, differentiation and embryonic development.
- Here we describe an *in vitro* fully kinetic neurite outgrowth assay miniaturised to 96 & 384-well microtitre plate formats based on analysis of time-lapse, phase-contrast images.
- Validation and pharmacology data from a range of cell types including human iPSC-derived neurons (iCell Neurons, CDI), primary neurons (rat cortex) and

- neuronal-like cell lines (Neuro-2A) are described. Neurite outgrowth (NOG) can be duplexed with cell count and cytotoxicity measurements using fluorescent probes.
- This approach affords a full temporal understanding of neurite outgrowth without the need for complex and expensive Ab-labelling methods (e.g. HCS).
- Miniaturisation of assays with iPSC-derived neurons (iCell Neurons, CDI) to 384-well format is a key step in maximising the value of these cellular reagents.



96 & 384-well NeuroTrack Assay – an integrated solution






Image processing & quantification




- Phase-contrast images and analysis masks show the kinetics of neurite outgrowth in rat cortical primary neurons.
- Time courses compare mean \pm SEM neurite length (mm/mm²) and branch point (1/mm²) values for rat cortical, human iCell Neurons and Neuro-2A cells (n=6, 24 and 48 wells respectively).

Validation: Cell body cluster normalisation




- Data represent mean \pm SD values for rat cortical and human iCell Neurons (n=6 wells).
- Note how normalization to cell body cluster allows for data comparison.

Nuclear-targeted GFP/RFP enables true cell counting




- Time courses for NuLight Green and Red Neuro-2A cells differentiated using atRA (18 μ M) showing the mean \pm SEM (A, B) (n=9 wells).
- Phase and Fluorescence images obtained for NuLight Green and Red Neuro-2A cells (C, D) 5 days post-treatment.

Kinetic multiplexed assays : NOG & neurotoxicity




NOG assay miniaturisation to 384-well format




- Cells treated with increasing concentrations of atRA (Neuro-2A) or BDNF (iCell Neurons) 1h post-plating and monitored for 5 and 4 days respectively. Increasing cell densities (iCell Neurons) in a 384-well plate format.
- Plate views show the kinetics of neurite length (mm/mm²) for Neuro-2A cells (A) or iCell Neurons (D).
- Inter-well variability plot; maximal neurite length \pm SEM (mm/mm²) for Neuro-2A (CV<7%, n=48 wells) (B).
- Bar graph shows the AUC of increasing iCell Neurons cell densities treated with increasing BDNF concentrations (CV=8%) (C) (n=24 wells).

Differential effects of Cytochalasin D across cell types



Quantitative pharmacology – Kinetic profiles



Acknowledgement

We thank Cellular Dynamics International for technical support with iCell Neurons.

